



Bed & Chair Exit Alarms



Our bed and chair exit alarms are designed to provide immediate alerts when a client leaves their bed or chair in order to alert carers to possible fall risks.

They can work as a stand alone system or are compatible with most nurse call systems with an optional connector lead.

The controller works in 1 of 3 selectable modes:

- 1 With a bed or chair sensor mat when a person attempts to get up from the bed or chair
- 2 With a floor sensor mat when a person steps on or applies weight to the sensor
- 3 With a magnetic pull cord that alerts staff if the client moves away from the controller e.g. attempts to leave a chair unsupervised



Works with bed, chair and floor sensor mats



Immediate or 2 second delay available



Magnetic clip alarms when client moves away




Silent, soft or loud audio levels can be selected



Silent mode available when connected to call systems



Battery powered with optional mains PSU



Standard Kit Contents:

- 1 × Fall Alarm Controller
- 1 × Bed Sensor Mat
- 1 × Wall Mounting Bracket
- 1 × Bed/Chair Mounting Bracket
- 1 × Magnetic Pull Cord
- 1 × 9V Battery

Optional Accessories:

- Chair Sensor Mat
- Heavy Duty Floor Sensor Mat
- Staff Pager
- Nurse Call Connection Lead



Distributed by:

Safety Systems Distribution Ltd
Unit 11B Hobson Industrial Estate
Burnopfield, Newcastle, NE16 6EA

Tel: 0800-3282950 Fax: 0800-3282951
www.safetysystemsdistribution.co.uk