

# Product Datasheet

## BC2

### Mini Breathcall Switch

- Pneumatic Breathcall Switch
- Just a puff of air to trigger
- Comes with protective saliva filter
- Adjustable clamp to fit side rail or headboard
- Flexible arm can be adjusted to fit any position
- Includes 10 Disposable filters
- Suitable for children and adults
- Links to directly to Intercall and most other Nursecall systems\*

\*Standard contact type: Normally open ¼" Mono Jack



### Features

#### Pneumatic

The Breathcall unit can be used as an aid to communication for persons with limited mobility, who are unable to use a standard call lead. For example: paraplegic and quadriplegic individuals and those with severe arthritic and rheumatic conditions.

#### Easy Positioning

The flexible tube means that the Breathcall unit can be positioned to suit the individual and can be easily transferred from one position to another.

#### Activation

Requiring only the lightest puff of air to activate, the Breathcall unit will then trigger the Alarm system or Nursecall system to alert the carer that assistance is required.

#### Operating Voltage

Up to 24 volts AC or DC at 33amps

#### Replacement Filters

10 clear soft vinyl straws are supplied.  
For additional Filters, use Order Code BCDF-10.



Please contact us for more details or to order:

Safety Systems Distribution Ltd

Unit 11 Hobson Industrial Estate, Burnopfield, Tyne & Wear, NE16 6EA

Tel: 0800-3282950 - Fax: 0800-3282951 - sales@nursecall.org - www.nursecall.org